

DIRECTIONS FOR USE: Correct Medication Labeling System

Pat. # US 6,955,002 B2; 7,381,288;
D538,851; D542,415

1 WRITE IN APPROPRIATE INFORMATION USING A PERMANENT MARKER

Afrin Strength: _____ % Expiration Date / Time: _____	Marcaine Strength: _____ % Expiration Date / Time: _____
Afrin Strength: _____ % Expiration Date / Time: _____	Marcaine Ephedrine Strength: <u>5%</u> / <u>200,000</u> Expiration Date / Time: <u>2-24</u> <u>1206</u>
Antibiotic 0.9% Saline Irrigation Name and Strength: _____ Expiration Date / Time: _____	Papaverine Strength: _____ % Expiration Date / Time: _____
Contrast Strength: _____ % Expiration Date / Time: _____	Contrast Strength: _____ % Expiration Date / Time: _____
Lidocaine Strength: <u>2%</u> Expiration Date / Time: <u>2-24</u> <u>1206</u>	Lidocaine Strength: <u>2%</u> Expiration Date / Time: <u>2-24</u> <u>1206</u>


2 ADHERE LABEL ONTO SYRINGE



3 PEEL PAPER BACKING HALFWAY BACK. ADHERE LABEL TO FLAG SURFACE. USE PRE-PRINTED OR MARK ON LABELS.

4 PEEL PAPER BACKING OFF TO EXPOSE ADHESIVE. COLORED SYRINGE LABELS WILL COME OFF WITH BACKING.


5 ADHERE FLAG TO OUTSIDE OF CUP. OPTIONAL: REMOVE ADHESIVE BACKING FROM BOTTOM OF CUP AND THEN ADHERE TO A SURFACE.




6 PEEL OFF SYRINGE LABEL. ADHERE CORRESPONDING COLORED LABEL TO SYRINGE BARREL.

7 SET SPECIMEN STRIP AND MARKING PEN ON BACK TABLE.

8 PLACE SPECIMEN ON SQUARE. WRITE SPECIMEN NAME DIRECTLY ON PAD. USING PERMANENT INK.




1 SPECIMEN ZONE	2 SPECIMEN ZONE	3 SPECIMEN ZONE	4 SPECIMEN ZONE
L _{LEFT} R _{RIGHT}	L _{LEFT} R _{RIGHT}	L _{LEFT} R _{RIGHT}	L _{LEFT} R _{RIGHT}



OPTIONAL: IDENTIFY LEFT OR RIGHT. USE NUMBERS OR LETTERS TO COUNT SPECIMENS.

REMEMBER TO TAKE A



Ansell Sandel Medical Solutions LLC
19736 Dearborn St. Chatsworth, CA 91311
818-534-2500 www.sandelmedical.com

Safety means NEVER having to say you're sorry.®

TIME OUT is a trademark of Sandel Medical Industries, LLC.

In the healthcare industry, the color orange is a trademark of Sandel Medical Industries, LLC.