


# DIRECTIONS FOR USE: DrifFloor™

Reorder No. 2419, 2440, 2449

# 1

PLACE DrifFloor™  
IN AREAS WHERE  
FLUID  
COLLECTS.  
WHITE SIDE  
FACES  
UP.



# 2

DISPOSE IN THE APPROPRIATE CONTAINER  
BY LIFTING THE MAT IN A MANNER THAT IT  
DOES NOT LEAK.

**Caution:** For safe use do not place on wet floor or in standing water. Do not over saturate to produce standing water. Use additional pads to absorb excess fluids.

( SINGLE USE ONLY )



Ansell Sandel Medical Solutions LLC  
19736 Dearborn Street, Chatsworth, CA 91311  
1-818-534-2500 [www.sandelmedical.com](http://www.sandelmedical.com)