

DIRECTIONS FOR USE:
TIME OUT Beacon™

Reorder No. 811-NS

**Take Time
For
TIME OUT®**


Use over the instrument tray, back table
or near the surgical site before incision.

The TIME OUT Beacon can be sterilized
as follows, according to your protocol:

- Gamma
- ETO (Ethylene Oxide)
- Steam

The TIME OUT Beacon is a
disposable product.

Do Not re-use.




SANDEL LLC
Leaders in Healthcare Safety®

Sandel Medical Industries, LLC
19736 Dearborn St. Chatsworth, CA 91311
818-534-2500 www.sandelmedical.com

Safety means NEVER having to say you're sorry.®

TIME OUT is a trademark of Sandel Medical Industries, LLC.

In the healthcare industry, the color orange is a trademark of Sandel Medical Industries, LLC.