


DIRECTIONS FOR USE: 2-in-1 Marker™

Reorder No. 1031 (Sterile)
Reorder No. 1031-NS (Non-Sterile)
Reorder No. 1031-L (Sterile)

Pat. # US D538,851 & Pat. # D542,415



Important!

1. Permanent Ink is not intended to be used on the skin.
2. For single patient use only.
3. Dispose of after use.
4. Skin surface must be clean and dried before marking the skin.
5. Alcohol or alcohol-based solutions will remove marks.

Safety means NEVER having to say you're sorry.®

TIME OUT is a trademark of Sandel Medical Industries, LLC.

In the healthcare industry, the color orange is a trademark of Sandel Medical Industries, LLC.


Leaders in Healthcare Safety®

Sandel Medical Industries, LLC
19736 Dearborn St. Chatsworth, CA 91311
818-534-2500 www.sandelmedical.com