

DIRECTIONS FOR USE: Single Tip Permanent Marker


Reorder No. 1000-L (Sterile)
Reorder No. 1000-N (Non-Sterile)

Pat. # US D538,851



- 1. Remove black cap for permanent marker before use.**
- 2. Replace cap after use.**

Permanent Ink.



Important!

1. Permanent Ink is not intended to be used on the skin.
2. For single patient use only.
3. Dispose of after use.


Leaders in Healthcare Safety®

Sandel Medical Industries, LLC
19736 Dearborn St. Chatsworth, CA 91311
818-534-2500 www.sandelmedical.com

Safety means NEVER having to say you're sorry.®

In the healthcare industry, the color orange is a trademark of Sandel Medical Industries, LLC.