



  
**SANDEL**  
Guard-It™

A simple alternative to the  
“scoop” technique!

#### ANTICIPATE THE UNEXPECTED

The newest addition to the SANDEL Engineered Sharps Injury Prevention line helps prevent accidental needlesticks by providing nurses and other healthcare professionals a way to safely cap a syringe needle.

This simple solution holds the syringe cap at an angle which makes capping the syringe with one hand easier and safer. This provides a simple alternative to the scoop method and complies with OSHA requirements\*.

*\*A feasible alternative to recapping*

#### PRODUCT INFORMATION

Item Number: 2280  
Packaging: 100/Cs (Sterile)  
Dimensions: 3" x 4" x 1/16"  
Disposable / Single Patient Use  
Compatible with standard syringes  
Orange Plastic (Recyclable\*)

*\*Must be disposed of per each facility's protocol.*

**Also available through your  
Custom Procedure Tray (CPT)  
Provider**

For more information, please visit:  
[www.sandelmedical.com](http://www.sandelmedical.com)  
or call 866-764-3327

