



Take Time For TIME OUT®

DIRECTIONS FOR USE: TIME OUT HOOD™

Patent Pending

Reorder No. 711

Place the TIME OUT Hood™ over the scalpel on the Mayo Stand or back table. The TIME OUT Hood™ is a reusable product.



CAUTION: Do not use alkaline-based detergents when cleaning this product as the solution will change the color of the aluminum

TIME OUT® is a registered trade mark of Sandel Medical Industries, LLC



Sandel Medical Industries, LLC
19736 Dearborn Street, Chatsworth, CA 91311 USA
1-818-534-2500 www.sandelmedical.com